

AUTOCLAVE LABELS

LABEL APPLICATIONS

Medical sterilization

Labels pre printed with a special color in order to monitor a correct sterilization process.

Available for:

- ❖ damp / steam sterilization
- ❖ ethylene oxide sterilization
- ❖ Gamma radiation on request


The indicator ink changes to a different color when the sterilization process is completed.

Documentation on the process can be provided by Checkpoint/Meto

Color change:


steam

pink -> brown




ethylene oxide

blue -> green



Gamma radiation

beige -> brown



The duplex label can be removed from the sterilized product and placed on the ID card of the patient. It is mainly used as a control label.

Sectors: institutions, pharmacies, laboratories, hospitals.

